

FFP3 NR  
Protective  
FOLDED MASK  
Q5M

The Purple's

16\*14\*15cm

The Purple's

FFP3 NR  
PROTECTIVE FOLDED MASK  
Q5M  
PFE > 99%



20 PCS

INDIVIDUALLY PACKED

EN149: 2001+A1:2009 FFP3 NR  
PPE-REGULATION(EU)2016/425

### Single-use face mask (Non-sterile)

#### Indications For Use

It is used for patients in order to reduce the risk of the spread of infection, and it is intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.

#### Limitations

- Not suitable for working places where open flames exist (such as welding, casting, etc).
- It is not suitable for professional protective use such as oxygen deficient, toxic gas environment.
- For healthcare use.

#### Checks prior to use

Please check whether the packaging is damaged, whether the mask body is damaged, deformed or has other obvious defects before use.

#### Fitting

1. Cup the mask in one hand.
2. Raise the mask to the face, ensuring that foam piece sits on the nose (if fitted). Use the other hand to ensure that the lower part is pulled down under the chin.
3. Hold mask firmly and place the lower strap over the head to the nape of the neck.
4. Place the upper strap on the crown of the head. Adjust straps to achieve a comfortable fit. Straps should not pass over the ears.
5. (if fitted) Form the nose clips around the nose. Hold the finger and thumb as shown, push the hand down so that the finger goes down one side of the nose and thumb the other. Maintain pressure against the nose to ensure that a seal is formed.
6. To check the fit of the mask, place both hands completely over the mask and exhale vigorously. If air leaks at the edge, readjust the mask, manipulating the nose bridge and straps to effect an acceptable seal.



If transporting this product use original packaging provided

The Purple's

FFP3 NR  
PROTECTIVE FOLDED MASK  
Q5M  
PFE > 99%



20 PC

INDIVIDUALLY PACKED

EN149: 2001+A1:2009 FFP3 NR  
PPE-REGULATION(EU)2016/425

#### Use cautions:

- Please check fit of mask prior to use;
- It is unlikely that the requirements for leakage will be achieved if facial hair passes under the face seal;
- Please do not use the mask under low air quality condition, such as serious contaminants in air, oxygen deficiency;
- Do not use the mask in explosive atmosphere;
- After use, please discard the mask following local regulation;
- The mask is only for single use, do not use repeatedly.


#### Storage Condition:

The product should be stored in a well ventilated, dark and dry environment. Keep away from fire, pollutants and possible pollution sources. Also, the transport storage 's temperature: -20 C ~38 C, relative humidity: <80%RH.

#### Shelf-Life:

The Protective Mask is valid for 3 years with the above storage condition.

Website: [www.thepurplesmask.com](http://www.thepurplesmask.com) sales@thepurplesmask.com  
Manufacturer: Zhejiang The Purples Protective Products Co.,Ltd.  
Address: No. 45, Suhua Street, Suxi Town, Yiwu, Zhejiang Province, China



Production Date: Refer to CQ  
Expired Date:  
Lot Number:  
6 973911 200415  
Made in China